
# **Iridium Catalyst for Hydrogenolysis of Urea Derivatives** to Formamides and Amines under Mild Conditions



(1,5-Cyclooctadiene)-[2-(2-diphenylphosphinophenyl)pyrrolido]iridium(l) 200ma

# Advantages

- Decomposes urea derivatives into formamides and amines
- Can be used under mild conditions
- Applicable on gram scale or with very low amount of catalyst (0.09 mol%)

# **Applications**

**Decomposition of symmetric urea derivatives** 

X = COOR, CONR2, CN, NR2, F, Cl, Br

C4021 can decompose urea derivatives into formamides and amines while retaining the reducing active functional group.

#### Decomposition of asymmetric urea derivatives

In the case of unsymmetrical urea derivatives, C4021 cleaves the C-N bond on the side with more substituents at nitrogen

# **Decomposition of polyurea**

C4021 shows highly catalytic activity against polyurea, and can decompose polyurea in a manner different from conventional

Reference T. Iwasaki, K. Tsuge, N. Naito, K. Nozaki, Nat. Commun. 2023, 14, 962. https://doi.org/10.1038/s41467-023-38997-2

#### **Related Product**

**Ru-MACHO®-BH** 

200mg / 1g [R0137]

For further information please refer to our website at www.TClchemicals.com. >>> TCI iridium





#### Ordering and **Customer Service**

#### TCI AMERICA

:800-423-8616 / 503-283-1681 Tel : 888-520-1075 / 503-283-1987 E-mail: Sales-US@TCIchemicals.com

### TCI EUROPE N.V.

: +32 (0)3 735 07 00 : +32 (0)3 735 07 01 E-mail: Sales-EU@TCIchemicals.com

#### **TCI Deutschland GmbH**

Tel : +49 (0)6196 64053-00 : +49 (0)6196 64053-01 E-mail: Sales-DE@TCIchemicals.com

#### Tokyo Chemical Industry UK Ltd. Tokyo Chemical Industry (India) Pvt. Ltd.

: +44 (0)1865 78 45 60 E-mail: Sales-UK@TCIchemicals.com

# 梯希爱(上海)化成工业发展有限公司

Tel :800-988-0390 / 021-67121386 : 021-6712-1385 E-mail: Sales-CN@TCIchemicals.com

#### TOKYO CHEMICAL INDUSTRY CO., LTD.

Tel: +81 (0)3-5640-8878

E-mail: globalbusiness@TCIchemicals.com

: 1800 425 7889 / 044-2262 0909

E-mail: Sales-IN@TCIchemicals.com

Condition A 4021 (3 mol%) H<sub>2</sub> (1 MPa) THF, 130°C

Condition B

H<sub>2</sub> (1 MPa)

Toluene, 130°C

C4021 (1-10 mol%) tBuOK (10-33 mol%)

• Chemicals itemized in this brochure are for research and testing use only. Please avoid use other than by chemically knowledgeable professionals. • Information such as listed products and its specifications and so on are subject to change without prior notice. • The contents may not be reproduced or duplicated in whole or in part without permission of Tokyo Chemical Industry Co., Ltd.