



HCS Scientific & Chemical Pte Ltd

www.hcs-lab.com

August Newsletter: Filtration

August Promotion: Vacuum Pump

The HCSVP050 and HCSVP050CR Vacuum Pump is suited for easy relocation due to its compact design. It can be used for various purposes such as filtration, gas transfers and more.



HCSVP050



HCSVP050CR
Corrosion Resistant

Terms & Conditions Applies*

Check out our latest listing of filters!

Our latest update are the Membrane filters and Syringe filters.


Click images to view the range we have to offer to you!



[Membrane filters](#)



Syringe filters [Simplepure](#)



Syringe filters [Chrompure](#)

For more information on the lab consumables that we have, you can visit [our site](#) or contact us with your list of product that you are looking for.

Is Syringe Filters new to you?

Syringe filters (sometimes called a **wheel filter** if it has a wheel-like shape) is a single-use filter cartridge. It is attached to the end of a syringe for use. Syringe filters may have Luer lock fittings, though not universally so. The use of a needle is optional; where desired it may be fitted to the end of the syringe filter.

A syringe filter generally consists of a plastic housing with a membrane which serves as a filter. The fluid to be purified may be cleaned by drawing it up the syringe through the filter, or by forcing the filtered fluid out, through the filter.

Syringe filters may be used to remove particles from a sample, prior to analysis done by HPLC. Particles can easily damage a HPLC system due to the narrow bore and high pressures within.

For more information on the [Simplepure and Chrompure](#), you can take a look at [the brochure](#).



Other **Filtration Apparatus** such as **Filter Papers** and **Filter Glasswares** are available at **HCS**

Follow our [social media](#) for future updates and information!



Follow us on our [social media](#) such as [Youtube](#), [Facebook](#), [Instagram](#), [Linkedin](#) to get more future informational updates such as, promotions, discounts, or events!

HCS SCIENTIFIC & CHEMICAL PTE LTD

Main Office: 2 Venture Drive #08-32 Vision Exchange S608526 | Warehouse: 2 Penjuru Place #01-08 2-8 Penjuru Tech Hub S608783



+65 6777 7077



+65 6777 0500



www.hcs-lab.com



marketing@hcs-lab.com

*You are receiving this email because you have an account with HCS Scientific & Chemical Pte Ltd.
Click here to unsubscribe to our mailing list.*