




HCS Scientific & Chemical Pte Ltd

Optical & Quartz Glassware

The Best Choice for your sampling



CUVETTES



RECTANGULAR CELLS

HCS Customizable Quartz Glassware Available!
For more information on the custom-made products Please contact us!

Other product range includes:

Crucible
Evaporating Dish
Volumetric Flask
Erlenmeyer Flask
Watch Glass
Beaker
Sight Glass
Quartz Pipe

Material	Transmission on Empty Cell	Deviation in Matching
Optical Glass	approx. 82% at 350nm	max 0.5% at 350nm
UV Quartz Glass	approx. 80% at 220nm	max 0.5% at 350nm
ES Quartz Glass	approx. 80% at 200nm	max 0.5% at 200nm
IR Quartz Glass	approx 88% at 2730nm	max 0.5% at 2740nm

HCS SCIENTIFIC & CHEMICAL PTE LTD

2 Penjuru Place #01-08, 2-8 Penjuru Tech Hub S608783

Tel: +65 6777 7077 Fax: +65 6777 0500 E-mail: marketing@hcs-lab.com

www.hcs-lab.com