


# SAFETY CABINETS

[ FOR FLAMMABLES ]



Standardized manufacture meets OSHA 29 CFR 1910 and NFPA 30 that is strict in accordance with EN14470-1 Approval, issued by German third-party authority. A number of features, such as patented door construction, heavy-duty steel, unique design of flush top, and UV-resistant epoxy coating etc, build up superior safety cabinets in the industry.




## 1 Quality epoxy coating

Excellent UV resistant, even covers the hidden areas



## 2 Anti-slippery handle

increase friction makes handling easy



## 3 Rounded safety corner

To reduce accidental cuts and potential injury.



## 4 Fully-welded construction

Reduce dirt collection and zero welding spot

Product	Capacity	Model	Style	Doors	Max Dimensions	Shelves	Weight
680411	4 / 15	Counter-top	Manual	1	560mm X 439mm X 438mm	1	26
690411	4 / 15	Counter-top	Self-close	1	560mm X 439mm X 438mm	1	28
681211	12 / 45	Compact	Manual	1	890mm X 599mm X 470mm	1	62
691211	12 / 45	Compact	Self-close	1	890mm X 599mm X 470mm	1	64
681511	15 / 57	Compact	Manual	1	1118mm X 599mm X 470mm	1	70
691511	15 / 57	Compact	Self-close	1	1118mm X 599mm X 470mm	1	71
681221	12 / 45	Piggyback	Manual	2	460mm X 1099mm X 470mm	0	59
691221	12 / 45	Piggyback	Self-close	2	460mm X 1099mm X 470mm	0	60
681721	17 / 64	Piggyback	Manual	2	610mm X 1099mm X 470mm	1	67
691721	17 / 64	Piggyback	Self-close	2	610mm X 1099mm X 470mm	1	68
6817211	17 / 64	Wall mount	Manual	2	610mm X 1099mm X 470mm	0	69
682021	20 / 76	Wall mount	Manual	2	1118mm X 1099mm X 305mm	5	93
682211	22 / 83	Stand alone	Manual	1	1660mm X 599mm X 470mm	3	98
692211	22 / 83	Stand alone	Self-close	1	1660mm X 599mm X 470mm	3	100
685411	54 / 204	Stand alone	Manual	1	1660mm X 599mm X 864mm	3	135
695411	54 / 204	Stand alone	Self-close	1	1660mm X 599mm X 864mm	3	136
682221	22 / 83	Undercounter	Manual	2	890mm X 899mm X 580mm	1	88
692221	22 / 83	Undercounter	Self-close	2	890mm X 899mm X 580mm	1	92
683021	30 / 114	Stand alone	Manual	2	1120mm X 1099mm X 470mm	1	103
693021	30 / 114	Stand alone	Self-close	2	1120mm X 1099mm X 470mm	1	107
683121	30 / 114	Stand alone	Manual	2	1660mm X 1099mm X 470mm	1	93
684521	45 / 170	Stand alone	Manual	2	1660mm X 1099mm X 470mm	2	148
694521	45 / 170	Stand alone	Self-close	2	1660mm X 870mm X 870mm	2	152
686021	60 / 227	Stand alone	Manual	2	1660mm X 870mm X 870mm	2	189
696021	60 / 227	Stand alone	Self-close	2	1660mm X 870mm X 870mm	2	193
689021	90 / 341	Stand alone	Manual	2	1660mm X 1099mm X 870mm	2	204
699021	90 / 341	Stand alone	Self-close	2	1660mm X 1099mm X 870mm	2	212



Social Medias

@hcslab



+65 6777 7077



+65 6777 0500



www.hcs-lab.com

**HCS Scientific & Chemical Pte Ltd**


2 Penjuru Place #01-08, 2-8 Penjuru Tech Hub, Singapore 608783

# SAFETY CABINETS

[ FOR CORROSIVES ]



The safety cabinets for corrosives, along with the safety cabinets for flammables, have the same standardized framework, same high quality material and sophisticated craftwork. The epoxy coating of the cabinets withstand corrosion from acids and bases. Moreover, the polypropylene trays on top of the adjustable galvanized shelves, collecting unexpected leakage and spill are corrosion resistant.



## 1 Quality epoxy coating

Excellent UV resistant, even covers the hidden areas



## 2 Polypropylene tray

Injected corrosion resistant tray is tenacious and durable



## 3 Rounded safety corner

To reduce accidental cuts and potential injury.



## 4 Fully-welded construction

Reduce dirt collection and zero welding spot

Product	Capacity	Model	Style	Doors	Max Dimensions	Shelves	Trays	Weight
680412	4 / 15	Counter-top	Manual	1	560mm X 439mm X 438mm	1	2	29
690412	4 / 15	Counter-top	Self-close	1	560mm X 439mm X 438mm	1	2	30
681025	10 / 38	Dual Storage	Manual	2	457mm X 1092mm X 457mm	0	1	73
681025	10 / 38	Dual Storage	Self-close	2	457mm X 1092mm X 457mm	0	1	74
681212	12 / 45	Compact	Manual	1	890mm X 599mm X 470mm	1	3	63
691212	12 / 45	Compact	Self-close	1	890mm X 599mm X 470mm	1	3	65
681512	15 / 57	Compact	Manual	1	1118mm X 599mm X 470mm	1	2	70
691512	15 / 57	Compact	Self-close	1	1118mm X 599mm X 470mm	1	2	71
681721	17 / 64	Piggyback	Manual	2	457mm X 1092mm X 457mm	1	1	67
691721	17 / 64	Piggyback	Self-close	2	457mm X 1092mm X 457mm	1	1	68
6817211	17 / 64	Piggyback	Manual	2	610mm X 1099mm X 470mm	0	0	69
682021	20 / 76	Piggyback	Self-close	2	1118mm X 1099mm X 305mm	5	5	93
682211	22 / 83	Stand alone	Manual	1	1660mm X 599mm X 470mm	3	3	98
692211	22 / 83	Stand alone	Self-close	1	1660mm X 599mm X 470mm	3	3	100
685411	54 / 204	Stand alone	Manual	1	1660mm X 599mm X 864mm	3	3	135
695411	54 / 204	Stand alone	Self-close	1	1660mm X 599mm X 864mm	3	3	136
682221	22 / 83	Undercounter	Manual	2	890mm X 899mm X 580mm	1	1	88
692221	22 / 83	Undercounter	Self-close	2	890mm X 899mm X 580mm	1	1	92
683021	30 / 114	Stand alone	Manual	2	1120mm X 1099mm X 470mm	1	1	103
693021	30 / 114	Stand alone	Self-close	2	1120mm X 1099mm X 470mm	1	1	107
683121	30 / 114	Stand alone	Manual	2	1660mm X 1099mm X 470mm	1	1	93
684521	45 / 170	Stand alone	Manual	2	1660mm X 1099mm X 470mm	2	2	148
694521	45 / 170	Stand alone	Self-close	2	1660mm X 870mm X 870mm	2	2	152
686021	60 / 227	Stand alone	Manual	2	1660mm X 870mm X 870mm	2	2	189
696021	60 / 227	Stand alone	Self-close	2	1660mm X 870mm X 870mm	2	2	193
689021	90 / 341	Stand alone	Manual	2	1660mm X 1099mm X 870mm	2	2	204
699021	90 / 341	Stand alone	Self-close	2	1660mm X 1099mm X 870mm	2	2	212



Social Medias

@hcslab



+65 6777 7077



+65 6777 0500



www.hcs-lab.com

**HCS Scientific & Chemical Pte Ltd**

2 Penjuru Place #01-08, 2-8 Penjuru Tech Hub, Singapore 608783

## SAFETY CABINETS

[ FOR HARSH CORROSIVES ]



The cabinets are built by high quality polypropylene. Which is harsh corrosion resistant. The compliant welding process makes easy cleaning and nice appearance, The sump at the bottom is 51mm in depth, collects unexpected leakage and spill.



Product	Capacity	Model	Style	Doors	Max Dimensions	Shelves	Weight
700415	4	Counter-top	Manual	1	559mm X 432mm X 432mm	0	13
701215	12	Stand alone	Manual	1	590mm X 460mm X 890mm	1	26
702215	22	Stand alone	Manual	1	590mm X 460mm X 1650mm	1	45
703025	30	Stand alone	Manual	2	1120mm X 1090mm X 460mm	1	65
704525	45	Stand alone	Manual	2	1660mm X 1090mm X 460mm	2	78

## SPECIALTY CABINETS

[ FOR TOXIC CHEMICALS ]



Most noxious gas is of higher density than air. High concentration will jeopardize lives. Given more accidents are due to the leakage of noxious gas, specialty cabinets deals with the toxic chemicals. Flexible electronic lock and accurate seamless design assure high level of safety.



Product	Capacity	Model	Style	Doors	Max Dimensions	Shelves	Trays	Weight
681214	12 / 45	Compact	Manual	1	890mm X 599mm X 478mm	1	2	63
682214	22 / 83	Stand alone	Manual	1	1660mm X 599mm X 492mm	3	4	99
683024	30 / 114	Stand alone	Manual	2	1120mm X 1090mm X 470mm	1	2	104
684524	45 / 170	Stand alone	Manual	2	1660mm X 1099mm X 495mm	3	4	149



## OIL WASTE CANS

[ FOR OIL WASTE DISPOSALS ]



A variety of features such as the innovating link rod of the lid, the unique noise reduction design, self-closing lid, and heavy-duty steel etc, makes the products robust and durable. The containers is always closed when not in use, so as to isolate oily disposals from fires and eliminate combustion risks.



Product	Capacity	Colour	Outer Diameter ( ØDiameter x Height)
280613	06 GL / 23 L	Red	Φ341mm X 383mm
280611	06 GL / 23 L	Yellow	Φ341mm X 383mm
281013	10 GL / 38 L	Red	Φ381mm X 471mm
281011	10 GL / 38 L	Yellow	Φ381mm X 471mm
281413	14 GL / 53 L	Red	Φ411mm X 552mm
281411	14 GL / 53 L	Yellow	Φ411mm X 552mm
282113	21 GL / 80 L	Red	Φ481mm X 579mm
282111	21 GL / 80 L	Yellow	Φ481mm X 579mm



## LIQUID PP CANS

[ FOR LIQUID WASTE DISPOSALS ]



These are made of high density PP material, which can resist strong corrosive chemicals. The tank body is equipped with self-close mechanism, preventing dents, penetration and fragility cause by corrosion. The can stops ignition source from entering the body while operating the waste can.

Product	Capacity	Outer Diameter ( ØDiameter x Height)
602515	2.5 L	Φ301mm X 370mm
605015	5 L	Φ301mm X 500mm

## SAFETY STORAGE CANS

[ FOR FLAMMABLE LIQUID ]



These cans provides security for flammable liquid storage and offer easy pouring to protect site personnel and minimize the risk of accident. The ergonomic design and production has taken safety into adequate consideration to meet up with the regulations.



Product	Capacity	Outer Diameter ( ØDiameter x Height)
582513	2.5 L	Φ294mm X 303mm
585013	5 L	Φ294mm X 446mm






## DRUM CABINETS

[ FOR FLAMMABLES ]



The drum cabinets are designed and built to the OSHA 29 CFR 1910.106 standard, and specialize in housing one or two drums. A particular feature stays with the double-drum roller. It is segregated into two units so that the drums does not affect each other during handling, and space utilization is also enhanced.



- 1 Quality epoxy coating**  
Excellent UV resistant, even covers the hidden areas



- 2 Optional ramp & roller**  
Easier unload and loading action




- 3 Shelf hanger**  
Offers "no-slip" stability, easy to adjust and leak proof



- 4 Fully-welded construction**  
Reduce dirt collection and zero welding spot

Product	Capacity	Style	Roller	Doors	Max Dimensions	Shelves	Weight
686031	1 - 55GL drums	Manual	No	2	1660mm X 870mm X 870mm	1	207
686131	1 - 55GL drums	Manual	Yes	2	1660mm X 870mm X 870mm	1	207
696131	1 - 55GL drums	Self-close	Yes	2	1660mm X 870mm X 870mm	1	211
689031	2 - 30GL drums	Manual	No	2	1660mm X 1099mm X 870mm	1	250
689131	2 - 30GL drums	Manual	Yes	2	1660mm X 1099mm X 870mm	1	250
699131	2 - 30GL drums	Self-close	Yes	2	1660mm X 1099mm X 870mm	1	254
689231	2 - 30GL drums	Manual	No	2	1660mm X 1686mm X 870mm	2	372
689331	2 - 55GL drums	Manual	Yes	2	1660mm X 1686mm X 870mm	2	372
699331	2 - 30GL drums	Self-close	Yes	2	1660mm X 1686mm X 870mm	2	376
689431	115GL, dual storage	Manual	Yes	2	1660mm X 1686mm X 870mm	2	372
699431	115GL, dual storage	Self-close	Yes	2	1660mm X 1686mm X 870mm	2	376



## DRUM CABINETS

[ ACCESSORIES ]

Product	Capacity	Max. Dimension	KG
67031	For 55GL x 1	749mm x 762mm	24
67032	For 30GL x 2	940mm x 762mm	25
67034	For 55GL x 2	1408mm x 762mm	45
67033	Ramp	711mm x 622mm	15



# SAFETY CABINETS

[ FOR COMBUSTIBLES ]




The safety cabinets for combustibles, along with the safety cabinets for flammables, have the same standardized framework, same high quality material and sophisticated craftwork. With more adjustable shelves the cabinets are adaptable to many common small containers for combustibles, such as spray cans and paint cans.



- 1 Quality epoxy coating**  
Excellent UV resistant, even covers the hidden areas



- 2 Anti-slippery handle**  
increase friction makes handling easy



- 3 Rounded safety corner**  
To reduce accidental cuts and potential injury.



- 4 Fully-welded construction**  
Reduce dirt collection and zero welding spot




- 5 Self-closing door**  
Mechanism concealed in the top layer provides obstruction-free access



- 6 Flame resistant vent**  
reduce risks of amassing noxious gas.



- 7 Safety handle**  
For extra security and safety



- 8 Adjustable leg**  
Provides stability on uneven surfaces

Product	Capacity	Model	Style	Doors	Max Dimensions	Shelves	Weight
681513	20 / 76	Stand alone	Manual	1	1180mm X 599mm X 470mm	2	71
691513	20 / 76	Stand alone	Self-close	1	1180mm X 599mm X 470mm	2	72
683023	40 / 151	Stand alone	Manual	2	1120mm X 1099mm X 470mm	3	108
693023	40 / 151	Stand alone	Self-close	2	1120mm X 1099mm X 470mm	3	113
684533	60 / 227	Stand alone	Manual	2	1660mm X 1099mm X 470mm	5	156
694533	60 / 227	Stand alone	Self-close	2	1660mm X 1099mm X 470mm	5	158
686023	120 / 454	Stand alone	Manual	2	1660mm X 870mm X 870mm	5	197
696023	120 / 454	Stand alone	Self-close	2	1660mm X 870mm X 870mm	5	199

## SAFETY CABINETS

[ FOR PESTICIDES ]



EN  
14470-1

Improper use of pesticides will cause great harm to workers. Filled with profession and care, presents the cabinets for this particular application. The cabinets contain solid steel walls and corrosion-resistant polyethylene trays. Not only effective for storage segregation of different pesticides, the cabinet also help to eliminate potential moisture and leakage.



Product	Capacity	Model	Style	Doors	Max Dimensions	PP Tray	Shelves	Weight
681216	12 / 45	Compact	Manual	1	890mm X 599mm X 470mm	2	1	63
691216	12 / 45	Compact	Self-close	1	890mm X 599mm X 470mm	2	1	65
683026	30 / 114	Stand alone	Manual	2	1120mm X 1099mm X 470mm	2	1	105
693026	30 / 114	Stand alone	Self-close	2	1120mm X 1099mm X 470mm	2	1	107
684526	45 / 170	Stand alone	Manual	2	1660mm X 1099mm X 470mm	3	2	148
694526	45 / 170	Stand alone	Self-close	2	1660mm X 1099mm X 470mm	3	2	150
686026	60 / 227	Stand alone	Manual	2	1660mm X 870mm X 870mm	3	2	189
696026	60 / 227	Stand alone	Self-close	2	1660mm X 870mm X 870mm	3	2	193
689026	90 / 341	Stand alone	Manual	2	1660mm X 1099mm X 870mm	3	2	204
699026	90 / 341	Stand alone	Self-close	2	1660mm X 1099mm X 870mm	3	2	212

## UNDER FUME HOOD CABINETS

[ FOR FLAMMABLES ]



EN  
14470-1

Based upon the solid experience in laboratory application, under fume hood cabinets with multiple models and functions. The toe-kick gives more utilization of undercounter space, and both sides are equipped with flame arrester vents.



Yellow	Grey	White	Capacity	Style	Doors	Max Dimensions	Shelves	Weight
681541	681544	681545	15/57	Manual, right hinge	1	908mm X 610mm X 549mm	1	63
691541	691544	691545	15/57	Self-close, right hinge	1	908mm X 610mm X 549mm	1	65
6815411	6815444	6815455	15/57	Manual, left hinge	1	908mm X 610mm X 549mm	1	105
6915411	6915444	6915455	15/57	Self-close, right hinge	1	908mm X 610mm X 549mm	1	107
681941	681944	681955	19/72	Manual	2	908mm X 762mm X 549mm	2	148
691941	691944	691955	19/72	Self-close	2	908mm X 762mm X 549mm	2	150
682341	682344	682345	23/87	Manual	2	908mm X 914mm X 549mm	2	189
692341	692344	692345	23/87	Self-close	2	908mm X 914mm X 549mm	2	193
683141	683144	683145	31/117	Manual	2	908mm X 1219mm X 549mm	2	204
693141	693144	683145	31/117	Self-close	2	908mm X 1219mm X 549mm	2	212



Social Medias

@hcslab



+65 6777 7077



+65 6777 0500



www.hcs-lab.com

**HCS Scientific & Chemical Pte Ltd**

2 Penjuru Place #01-08, 2-8 Penjuru Tech Hub, Singapore 608783

## PPE CABINETS

[ FOR PROTECTIVE EQUIPMENTS ]



**EN**  
14470-1

Designed with two compartments for chemical protective clothing hung on aluminium alloy rails, and other personal protective equipments



**1**

### Cabinet hinges


for safety with good corrosion resistant



**2**

### Safety hammer

Easy to acquire without any hamper



Product	Capacity	Style	Doors	Max Dimensions	Roller	Shelves	Weight
984521	45	Manual	2	1800mm X 1000mm X 600mm	No	3	87

## STORAGE CABINETS

[ FOR REAGENTS ]



**EN**  
14470-1

made of high grade PP boards of 8mm, corrosion resistant and leakage-proof resulting from seamless welding process. Internal storage compartment is adjustable. Tempered glass window and double lock design



Product	Capacity	Style	Doors	Max Dimensions	Roller	Shelves	Weight
714025	40L	Manual	2	1800mm X 900mm X 450mm	No	4	60
714045	40L	Manual	4	1800mm X 900mm X 450mm	No	2	60
7140245	40L	Manual	4	1800mm X 900mm X 450mm	No	2	60



Social Medias

@hcslab

**HCS Scientific & Chemical Pte Ltd**

2 Penjuru Place #01-08, 2-8 Penjuru Tech Hub, Singapore 608783



+65 6777 7077



+65 6777 0500



www.hcs-lab.com